



# Self-collection Guide for HPV DNA Test




1 Gather swab, bottle with clear solution, and biohazard bag from the kit




2 Wash your hands with soap




3 Turn cap to remove swab




4 Hold swab at the red marking, do not touch the tip of swab



5 Position yourself for sample collection



6 Insert swab until fingers are at vagina opening



7 Rotate swab 10 to 30 seconds



8 Uncap bottle with clear solution



9 Swirl swab with it fully immersed 20 seconds



10 Drain fluid from swab



11 Close the bottle tightly, place it in the biohazard bag, and discard swab



12 Store the bottle in room temperature and away from direct sunlight



13 Return bottle to clinic within 72 hours